



## SPEEDBLOCKS HEAD IMMOBILIZER

### TECHNICAL CHARACTERISTICS

#### PACKAGING

- 1 BASE
- SIDE RISERS
- STRAP AND BEARING SET

#### USE AND STORAGE TEMPERATURE RANGES:

WHEN THE STOP IS PUT AWAY: FROM -34 TO 52 .

WHEN USING THE STOPPER: FROM -18 TO 43

#### CODICE PRODOTTO

210000440 SPEEDBLOCKS HEAD IMMOBILIZER

210000441 BLOCKS FOR SPEEDBLOCK HEAD STOPS

210000442 STRAP SET AND BEARINGS FOR SPEEDBLOCK ( 5 PZ)

983060 BASE FOR SPEEDBLOCK HEAD IMMOBILIZER LAERDAL



Speedblocks Starter Pack for BaXstrap Backboard. The SpeedBlocks offers a double advantage: despite being designed for durability and performance, it has the price of a disposable item.

The application of the SpeedBlocks is carried out with a simple four-step procedure: an inimitable multi axial locking system allows adaptation to every shape of head and position of wounds as well as direct access to the patient's ears for a better inspection and communication.

The inexpensive blocks with disposable foam pads and head straps are so effective that you no longer have to secure the patient's head with adhesive tape.

#### Attention:

The head immobilizer SpeedBlocks should be used on patients from 2 years of age up and is not suitable for vacuum mats.

SpeedBlocks should only be used by specially trained people.

#### TECNOLIFE S.r.l. -

Via Cimarosa, 71 - 40033 Casalecchio di Reno (BO)

tel.051/6259059 – fax.051/0510419

info@tecnolifesrl.com - www.tecnolifesrl.com - www.tecnolifeshop.com