



LAERDAL THOMAS TUBE HOLDER

TECHNICAL CHARACTERISTICS



NON STERILE PRODUCT



SINGLE USE PRODUCT



LATEX FREE

USE TIME MUST NOT EXCEED 24 HOURS

OPERATING TEMPERATURE: - 18 °C UP TO 43 °C

TEMPERATURE FOR STORAGE:
- 34 °C UP TO 52 °C

PACKAGING

PACK OF 1 PC

PACKAGE OF THE THOMAS TUBEHOLDER

FOR ADULTS:

PACK OF 25 PCS

CERTIFICATIONS

CND: R0199

RDM: 109229

PRODUCT CODE

210000454 LAERDAL THOMAS TUBE
HOLDER ADULT MODEL

210000454. LAERDAL THOMAS TUBE
HOLDER SELECT ADULT MODEL
(PACK OF 25 PCS)



210000454 ADULT



210000454. ADULT SELECT

Designed to provide a range of advanced respiratory devices, the Thomas tube holders are easy to apply and can provide a reliable grip compared to tape fixing methods. The unique and improved design of the Thomas tube holder makes it possible to safely stop endotracheal tubes, laryngeal masks, combitubes, LTA and other supraglottis airway devices.

The device is equipped with a Quick-set screw system that keeps the pipes securely positioned and facilitates the adjustment or removal.

Thanks to a large opening it allows the application of a wide range of single or double way tubes.

It allows a fast and efficient application of the strap around the patient's neck without head movement.

The large opening also allows easy suction and better inspection of the mouth. The velcro strap provides safety and comfort.

Depending on the color code, the adult and pediatric devices can be more easily identified.

Ventilations can continue throughout the application.

The pediatric and adult model allow the accommodation of external diameters between 4.3 mm to 21 mm.

Concerning the Thomas Select tube holder device, it is equipped with a wider opening than the Thomas tube holder. It allows to fix a more numerous range of supraglottic devices for the respiratory tract, laryngeal masks and endotracheal tubes. Thanks to the push screws, the tubes can be fixed quickly and efficiently.

TECNOLIFE S.r.l. -

Via Cimarosa, 71 - 40033 Casalecchio di Reno (BO)

tel.051/6259059 - fax.051/0510419

info@tecnolifesrl.com - www.tecnolifesrl.com

Comparison of the two models: Thomas Tube holder Adult model and Thomas Tube holder Select Adult model

	Thomas Tube Holder	Thomas Select Tube Holder
		
Tube housings measurements	≥6.5 mm(internal diameter) up to 21 mm (external diameter)	≥ 6.5 mm (internal diameter) up to 42 mm (external diameter)
Mechanism with push screws : for a quick initial adjustment on the tube locking device		✓
Screw mechanism to secure it to the pipe locking device	✓	✓
Release button for quick release		✓
Wide opening	✓	✓
Nozzle for the suction	✓	✓
Wide band for closure around the head	✓	✓
Flexible anti-bear device	✓	✓
Symbol for the right orientation on the nose, is located above the opening	↔	✓

TECNOLIFE S.r.l. -

Via Cimarosa, 71 - 40033 Casalecchio di Reno (BO)

tel.051/6259059 – fax.051/0510419

info@tecnolifesrl.com - www.tecnolifesrl.com